




Search



WE ARE A LEADING MANUFACTURER OF MICROSCOPE

Products

- Biological Microscope
- Stereo Microscope
- ZOOM
- XTL**
- XTL6445
- XTL6745
- HSM
- GEM
- DSM
- Accessory
- Metallurgical Microscope
- Polarization Microscope
- Fluorescence Microscope
- Industry Microscope
- Image Analysis & Camera
- Accessory

XTL6445



**XTL6445--- Streamline Integrated Mechanical Design**

Easy operation, long working distance, clear resolved image, enclosure integrated.

It is used extensively world medical research and health care, biology ,botany research and agriculture, as well as in electronic component manufacturing.

It is specially suited for LED, PCB inspection, pressing plating and electronic component inspection.

The XTL6445 with 0.7x~4.5x zoom magnification, inclination angles of 45° respectively.

Specification	XTL6445
<b>Microscope Body-Standard</b>	
Optical System	Parallel Optical System
Head Observation	Binocular 45°
Interpupillary Adjustment	55~75mm
Zoom ratio	1:6.4
Magnification Range	0.7X~4.5X
Objective	1X
Working Distance	100mm
Eye-piece	WF10X/F.N.20
Field of View	28.6~4.4mm
<b>Microscope Body-Optional</b>	
Magnification Range	3.5X~180X
Auxiliary Objectives	0.5X, 1.5X, 2X
Working Distance	30~165mm
Eye-piece	WF15X/F.N.15, WF20X/F.N.10
Field of View	1.1~57.1mm
<b>Microscope Accessories</b>	
Illumination	Fluorescent Right light, Cold Light Source Series
Digital Cameral	DCM, MDC series
Adapter (Trinocular)	CCD Adapter Series

XTL6445B-B1 XTL6445T-B1	XTL6445B-B2 XTL6445T-B2	XTL6445B-B3 XTL6445T-B3	XTL6445B-B4 XTL6445T-B4	XTL6445B-B5 XTL6445T-B5
XTL6445B-J1 XTL6445T-J1	XTL6445B-J2 XTL6445T-J2	XTL6445B-J3 XTL6445T-J3	XTL6445B-J4 XTL6445T-J4	
XTL6445B-WZ1 XTL6445T-WZ1	XTL6445B-WZ2 XTL6445T-WZ2			