



**MetaDi® Supreme Polycrystalline Diamond Pastes**

**MetaDi® Supreme Polycrystalline Diamond Pastes** offer the same benefits as our MetaDi® Supreme Suspensions in paste form. Polycrystalline diamonds have a greater number of cutting facets than the monocrystalline diamonds found in MetaDi® or MetaDi® II Diamond Paste. The increase in cutting facets provides higher removal rates and reduced subsurface deformation. This product is naturally colored with no dyes.

### MetaDi® Supreme Polycrystalline Diamond Paste, Heavy Concentration Diamond

Diamond Size	Concentration	Tip Color	Quantity	Number	Price
1 micron	Heavy	Blue	5g	40-6300	
1 micron	Heavy	Blue	20g	40-6301	
3 micron	Heavy	Green	5g	40-6302	
3 micron	Heavy	Green	20g	40-6303	
6 micron	Heavy	Yellow	5g	40-6304	
6 micron	Heavy	Yellow	20g	40-6305	
9 micron	Heavy	Deep Red	5g	40-6306	
9 micron	Heavy	Deep Red	20g	40-6307	
15 micron	Heavy	Brown	5g	40-6308	
15 micron	Heavy	Brown	20g	40-6309	



**Replicating Tape, Replicating Adhesive Bases, Disc Holders, and Flapper Wheel**

**Replicating Materials** are primarily used for field metallographic applications and subsequent laboratory analysis.

Description	Quantity	Catalog Number	Price
Disc Holder, w/ groove, hard rubber 1" (25.4mm)	1	57-4550	
Disc Holder, w/ groove, sponge rubber 1" (25.4mm)	1	57-4552	
Disc Holder, hard rubber 1" (25.4mm)	1	57-5005	
Disc Holder, hard rubber 2 7/8" (73mm)	1	57-5003	
Replicating Tape 1" x 0.005" x 25' (25.4mm x 0.13mm x 7.6m)	1 roll	57-5000	
Replicating Adhesive Bases 1" x 3" (25.4mm x 76.2mm)	30 per box	57-5001	
Replicating Kit*	1	57-5002	
Flapper Wheel, 60 (P60) grit	1	57-5004	
Manual Applicator for Replicating Media	1	57-6000	
Static Mixers	10	57-6001	
Replicating Media 1/7, 50ml (1.7oz) 1 minute working time, 7 minute cure time	1	57-6002	
Replicating Media 7/30T, 50ml (1.7oz) 7 minutes working time, 30 minute cure time	1	57-6003	
Replicating Set**	1	57-6100	



**Replicating Media**

\* Replicating kit includes 10 each 3" x 3" (76.2mm x 76.2mm) Replicating Sheets and Adhesive Bases

\*\* Replicating set includes 1 each of catalog numbers 57-6000, 57-6001, and 57-6002