




Search



WE ARE A LEADING MANUFACTURER OF MICROSCOPE

Products

Biological Microscope

Stereo Microscope

ZOOM

ZOOM460

ZOOM645

ZOOM1020

ZOOM630

Accessory

XTL

HSM

GEM

DSM

Accessory

Metallurgical Microscope

Polarization Microscope

Fluorescence Microscope

Industry Microscope

Image Analysis & Camera

Accessory

ZOOM645



ZOOM645--- Airtight, Anti-mold and Anti-electrostatic stereomicroscope.

This perfect microscope was developed specially to meet the needs of customers who demand more advanced features and higher performance at lower prices.

The Airtight, Anti-mold and Anti-electrostatic ensures peace of mind under any operating environment.

The ZOOM645 with 0.8x~5x zoom magnification, inclination angles of 45° respectively. It can be with 0.5x and 2x auxiliary lens.

Specification	ZOOM645
Microscope Body-Standard	
Optical System	Greenough stereoscopic
Head Observation	Binocular, 45° inclined
	Trinocular 45°, inclined with C adapter
Interpupillary Adjustment	52~76mm
Diopter Adjustment	±5°
Zoom ratio	Zoom645 1:6.3
Magnification Range	0.8X~5X
Objective	1X
Working Distance	115mm
Eyepiece	WF10X/F.N.22
Field of View	27.5~4.4mm
Microscope Body—Optional	
Magnification Range	4X~200X
Auxiliary Objectives	0.5X, WD = 211mm
Auxiliary Objectives	2X, WD = 43.5mm
Working Distance	43.5~211mm
Eyepiece	WF15X/F.N.16, WF20X/F.N.12.5
Field of View	2.2~55mm
Accessory—Optional	
Illumination	Fluorescent right light
Digital Cameral	DCM, MDC series
C Adapter(Trinocular)	CDD Adapter Series

Pack:

Only Head: 1 carton, G.W.: 2.5KGS, Meas.: 25*25*24CM

Head with Stand: 1 carton G.W.: 9.5KGS (with LED) 8.5KGS (without LED), Meas.: 36*30*40CM