

BUEHLER® Delta® PetroCut™ GEOLOGICAL CUTTER

- Affordable, Rugged Manual Cutter Designed for Rocks, Minerals, Nonmetallic and Geological Specimens
- Features Universal Clamping with a Rock Clamp Assembly
- Produces Serial Cuts Quickly and Consistently
- X-Y Table Bed
- Specimen Sizes Up To 3.75" (95mm) in Diameter
- Uses 10" (254mm) Diamond Blade
- Compact Tabletop Design
- Heavy-Duty Aluminum Base
- Large Viewing Window and Lighted Cutting Chamber
- Clean-Out Hose



BUEHLER



The Delta® PetroCut™ Geological Cutter is a manual, 4Hp (3.0kW) cutter using a 10" (254mm) Continuous Rim Diamond Blade or Rimlock* Diamond Blade for sectioning of rock, minerals, ceramics, coal, glass, concrete, refractory and many geological samples up to 3.75" (95mm) in size. Features include a Rock Clamp Assembly with an X-Y table bed and cross-feed table for serial cutting. Clamps can be positioned at various points on the table bed. Other specimen vises are available for holding irregular shapes.

The Delta PetroCut Geological Cutter includes a clean-out hose for quick debris removal. Two jets optimize presentation of the coolant to the wheel and specimen. Safety is maintained with an emergency stop, hood lock and a mechanical brake which is engaged when the power is turned off. The heavy duty aluminum cast base provides stability and corrosion resistance. Operation instructions are included.

10-2156-260 Delta® PetroCut™, 230 VAC, 60Hz, 3 Phase

10-2156-460 Delta® PetroCut™, 460 VAC, 60Hz, 3 Phase

10-2156-250 Delta® PetroCut™, 220 VAC, 50Hz, 3 Phase

10-2156-400 Delta® PetroCut™, 380 VAC, 50Hz, 3 Phase

Technical Data

Operation: Manual Chop Cutting

Motor Power: 4.0Hp (3.0kW)

Blade Size and Speed: 10" (254mm) Diameter Blades, 3400 rpm (60Hz), 2830 rpm (50Hz)

Cutting Envelope: Sample Sizes up to 3.75" (95mm)

Controls: Push Button Start and Stop/Emergency Stop; Coolant automatically activates when cutter is started

Other Features: Rock Clamp Assembly, Clean-out Hose, Heavy Duty Interior Lights (14W), Large Viewing Window

Shipping Weight: 200 lbs. (91kg)

Dimensions: 24" W x 17" H x 26.5" D (600mm x 425mm x 665mm)

Accessories

Blade Options **

11-4710 Continuous Rim Diamond Blade (for geological materials - yields best quality cuts) 10" x 0.060" x 1¼" (254mm x 1.5mm x 32mm)

11-4120 Rimlock* Diamond Blade (for geological materials - yields good quality cuts) 10" x 0.040" x 1¼" (254mm x 1.016mm x 32mm)

Optional Accessories

10-2165-260 Recirculating System, 230 VAC, 60Hz, 3 Phase

10-2165-460 Recirculating System, 460 VAC, 60Hz, 3 Phase

10-2165-250 Recirculating System, 220 VAC, 50Hz, 3 Phase

10-2165-400 Recirculating System, 380 VAC, 50Hz, 3 Phase

16-1497 Recirculating System Cart, with Wheels

Vise Accessories

10-3525 Irregular Specimen Speed Vise Kit (Large)

10-3531 Irregular Specimen Speed Vise Kit (Small)

10-3533 Rock Clamp Assembly



Delta® PetroCut™ Cutter showing Rock Clamp Assembly and Clean-out Hose.

Rock Clamp Assembly with specimen.

Optional Speed Vise Accessory Kit allows increased flexibility when clamping irregular shapes.

*Rimlock is a registered trademark of Felker Operations **Blade not included with Delta® PetroCut™

For a complete listing of Buehler Equipment and Consumables, please refer to our Buehler Equipment and Consumables Buyer's Guides. Buehler continuously makes product improvements; therefore, technical specifications are subject to change without notice.

© 2006 BUEHLER LTD. Printed in U.S.A. 25M0106 FN01289

*For metallurgical equipment produced by BUEHLER LTD. USA

BUEHLER ANALYST® SECTION

6



BUEHLER

BUEHLER LTD. - Worldwide Headquarters
41 Waukegan Road • P.O. Box 1
Lake Bluff, Illinois 60044-1699 USA
Tel: 847/295-6500 • Fax: 847/295-7979
Sales: 1/800/BUEHLER • 1/800/283-4537
Web Site: <http://www.buehler.com>
Email: info@buehler.com

BUEHLER GMBH

In der Steele 2 • 40599 Düsseldorf
Postfach 16 03 55 • 40566 Düsseldorf
Telefon: (+49) 0211/974100 • Telefax: (+49) 0211/9741079
Web Site: <http://www.buehler-met.de>
Email: info@buehler-met.de

BUEHLER SARL

Miniparc de Dardilly, Bât. 0
3, Chemin du Jubin • 69570 Dardilly-France
Téléphone: (+33) (04) 37 59 81 20
Télécopie: (+33) (04) 37 59 81 29
Web Site: <http://www.buehler.fr>
Email: info@buehler.fr

BUEHLER UK

101, Lockhurst Lane • Coventry • CV6 5SF
Telephone: (+44) (0) 24 7658 2158 • Fax: (+44) (0) 24 7658 2159
Web Site: <http://www.buehler.co.uk>
Email: sales@buehler.co.uk

BUEHLER CANADA

9999 Highway 48
Markham, Ontario L3P 3J3
Telephone: (905) 201-4686 • Fax: (905) 201-4683
Sales Telephone: 1-800-268-3593
Email: info@buehler.ca

BUEHLER ASIA

5/F Vogue Centre
696 Castle Peak Road
Lai Chi Kok, Kowloon
Hong Kong, SAR, China
Telephone: (852) 2307 0909
Fax: (852) 2307 0233