

BUEHLER® DELTA™ VISE KITS



10-2170 Sliding Vise Left and 10-2171 Sliding Vise Right are versatile quick release style vises for routine clamping of large to small specimens. Used with 10-2100 series DELTA Abrasive Cutters and 95-C1800 series ABRASIMET® Abrasive Cutters.



95-C1821 MetKlamp VIII Cam-Lock Vise Left and 95-C1822 MetKlamp VIII Cam-Lock Vise Right for routine clamping of medium to small specimens. Used with 10-2100 series DELTA Abrasive Cutters and 95-C1800 series ABRASIMET Abrasive Cutters.



10-2245 MetKlamp VII Vise Left and 10-2246 MetKlamp VII Vise Right are versatile cam-lock style vises for routine clamping of large to small specimens. Used with 10-2300 series DELTA Abrasive Cutters.



10-3523 Irregular Specimen Speed Vise Kit provides a vertical clamping force for large height and irregular shaped specimens. Used with 10-2200 series and 10-2300 series DELTA Abrasive Cutters.



10-3526 Horizontal Clamp Speed Vise Kit allows the Speed Vise Kit (10-3523) to be used horizontally. This is particularly helpful in clamping parts requiring longitudinal sectioning. Used with 10-2200 series and 10-2300 series DELTA Abrasive Cutters.



10-3525 Adjustable Vee Blocks provide a variable vee location and support for round and irregular shaped specimens. Used with 10-2200 series and 10-2300 series DELTA Abrasive Cutters.



10-3527 Vee Block Clamp Kit cradles specimens in a fixed vee. Stainless steel shims are included to assist in height adjustment. Used with 10-2200 series and 10-2300 series DELTA Abrasive Cutters.



10-3528 Large Riser Block Assembly Kit is used under the 10-3523 Speed Vises to increase their height for clamping of large specimens. **10-3532 Medium Riser Block Assembly Kit (not shown)** is used under the 10-3531 Speed Vises to increase their height for clamping of medium specimens. Both are used with 10-2200 series and 10-2300 series DELTA Abrasive Cutters.



10-3529 Sliding Vise Left and **10-3530 Sliding Vise Right** are versatile quick release style vises for routine clamping of large to small specimens. Used with 10-2200 series and 10-2300 series DELTA Abrasive Cutters.



10-3531 Irregular Specimen Speed Vise Kit provides a vertical clamping force for medium height and irregular shaped specimens. Used with 10-2200 series and 10-2300 series DELTA Abrasive Cutters.

Buehler continuously makes product improvements; therefore, technical specifications are subject to change without notice.

© 2002 BUEHLER LTD. Printed in U.S.A. 20M1002 FN01198 *For metallurgical equipment produced by BUEHLER LTD. USA

BUEHLER ANALYST[®]
SECTION

12



BUEHLER

BUEHLER LTD. – *Worldwide Headquarters*
41 Waukegan Road • P.O. Box 1
Lake Bluff, Illinois 60044-1699 USA
Tel: 847/295-6500 • Fax: 847/295-7979
Sales: 1/800/BUEHLER • 1/800/283-4537
Website: <http://www.buehler.com>
Email: info@buehler.com

BUEHLER GMBH
In der Steele 2 • 40599 Düsseldorf
Postfach 16 03 55 • 40566 Düsseldorf
Telefon: (+49) 0211/974100 • Telefax: (+49) 0211/9741079
Web Site: <http://www.buehler-met.de>

BUEHLER SARL
68, Chemin des Ormeaux
F-69760 Limonest
Téléphone: (+33) (04) 72 17 92 50 • Télécopie: (+33) (04) 72 17 92 49

BUEHLER UK
Milburn Hill Road • University of Warwick Science Park
Coventry CV4 7HS United Kingdom
Tel: (+44) (0) 24 7669 2242 • Fax: (+44) (0) 24 7669 3032
metsales@buehler.co.uk